
<b>de</b>	Gebrauchsanweisung: Stimmgabeln.....	2
<b>en</b>	User's Manual: Tuning Forks .....	5
<b>fr</b>	Mode d'emploi : Diapasons .....	8
<b>it</b>	Istruzioni per l'uso: Diapason .....	11
<b>es</b>	Instrucciones de empleo: Diapasones.....	14
<b>pt</b>	Manual de operação: Diapasões .....	17
<b>ru</b>	Руководство по применению: Камертоны ..	20

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
## User's Manual Rydel Seiffer, Lucae and Hartmann Tuning Forks

Example of Rydel Seiffer tuning fork,  
c 128/C 64 Hz with damper and  
plastic foot

Dear Customer, thank you for choosing a KaWe product. Our products are known for their high quality and longevity. This KaWe product meets EC Standards 93/42/EWG (standards for medical products).

**Please read this User's Manual thoroughly and carefully before attempting to use this product and heed the given care instructions.**

**Use:** These tuning forks are only to be used by authorized, trained personnel. The tuning forks create sound that is then conducted through the air and bone. The tuning forks are only to be placed on healthy skin. The plastic foot or end of the vibrating fork is placed on to the part of the body that is being examined. These are primarily spots on the extremities where the joints are not covered by muscle or fatty tissue. In order to ensure long-term monitoring, the results of each inspected area should each be compared with one another over time. For this purpose, the response of the patient is needed in order to determine at which point of the fading process of the fork's oscillation the patient no longer feels the vibration.

**Using the Rydel Seiffer tuning fork:** The Rydel Seiffer tuning fork produces two different frequencies: with attached dampers, the tone C = 64 Hz (as printed on the dampers) and the tone C = 128 Hz (as printed on the tuning fork itself). If the dampers are attached, the bottom edge of the dampers must be set to the calibration mark on the tuning fork tines whereby the label "c 64" on the dampers and the label "c 128" on the tuning fork (front side) as well as the calibration mark are to always be facing the observer. The dampers can be loosened and/or tightened with the knurled screws. Please tighten these knurled screws securely.

There are scales on the dampers on the back side of the tuning fork. By observing these scales, the vibration intensity of the tuning fork can be started with the same magnitude of vibration each time the patient is examined.

This takes place as follows:

- The tuning fork with attached dampers (thus/so the tone C = 64 Hz) is struck on the ball of the hand. The initial vibration of the tuning fork tines is so strong that both triangles appear to be completely blurred. After a few seconds, depending on the intensity with which the tuning fork was struck, the triangles will begin to appear on the base line and will grow slowly, that is with slowly diminishing vibration intensity to the first, second, third etc. tick mark on the scale. This means that the examination can be begun on a particular tick mark on the scale, whereby, depending on the light quality and/or visibility, the black or white triangle can be used.

It is well-known that the onset of the deterioration of the nervous system can be detected through a decrease in vibration perception. A replication of the magnitude of these vibrations is a prerequisite for the practical application of this diagnosis.

**Using the Lucae tuning fork:** The Lucae tuning fork produces 7 different frequencies from c 128 to b 240. The exact frequencies are given in the following table.

	c	d	e	f	g	a	b
Frequencies in Hz	128	144	160	170,66	192	213,33	240

The frequency is changed by sliding the weights on the tines of the fork. The weights can be easily moved with the hand. When changing the tone it is absolutely imperative that the spring in the weight slides securely into the slot on the tuning fork. Both weights must always be set at the same height. It is not possible to create overtones by setting

one weight to a different tone than the second one. The tone produced is designated by the stamped letter, which is found directly below the weight.

The tuning fork is to be struck preferably with the heel of the hand or with a piece of wood or plastic, but never with an object that is too hard (such as metal, steel or stone). The tuning fork should be struck in a brief and spring-like manner, meaning the tuning fork is only touched very briefly. Striking the tuning fork in this brief manner ensures the longest possible duration of vibration.

During the examination (sensitivity test etc.) the tuning fork is struck and then the base of the black plastic foot on the tuning fork is placed on the part of the patient's body that is to be examined.

**Purpose:** During hearing, airborne sound normally travels through the outer ear channel to the eardrum. The signal is transferred by the ossicle of the middle ear to the inner ear and is there taken in by the receptor cells (aerotympanal conduction). Sound oscillations can also be transferred through the bones, by-passing the middle ear, directly to the inner ear (bone conduction). The transfer of sound oscillations through air channels or bone channels depends on the frequency and differs between aerotympanal conduction and bone conduction.

**Unsuited use/contraindication of the product:** Any use other than that described here is not in accordance with the intended use of the unit. The manufacturer is not liable for any resulting damages. The user alone bears the risk.

**Further information, maintenance, storage:** When used and stored properly, this product will serve you many years.


**Guarantee:** When used properly and in accordance with the User's Manual we guarantee this product for two years from the date of purchase. Should you need further information or should your instrument require repair, please contact your dealer.

**Warning:** Please keep in mind that even minor damage to the tuning fork can influence the oscillation frequency. If the tuning fork has weights that are securely screwed on to the tines, these may not be unscrewed. The frequency will be changed if the weights are moved and the exact positioning of the weights can only be carried out by the manufacturer.

**Limitation of reuse preparation:** Frequent preparation has only a minimal effect on the instrument. The end of the product's life is normally due to wear and tear and damage by use. It is recommended that these instruments be re-prepared as soon as possible after their use.

## Instructions

**Storage and transport:** Rydel Seiffer tuning fork: adjusting the tuning fork after transport: During transport it is possible for the dampers to become loose and move. Before using the tuning fork ensure that both dampers are fastened on the tuning fork such that the markings on the back side of the tuning fork line up with the lower edge of the damper (see picture). Only then can we guarantee your proper functionality (frequency) of the tuning fork.

**Cleaning preparation:** No special requirements.

**Manual cleaning:** After use, the tuning forks may not be cleaned in a steriliser. To clean the tuning fork and the foot, gently wipe both with a dry or slightly greased cloth.

**Automatic cleaning:** Please ensure that the cleaning system is suitable for medical instruments. Only cleaning agents suitable for metals are to be used. The concentrations and reaction times given by the manufacturer of the cleaning agents are to be followed.

**Disinfection:** Parts that come into contact with the patient are to be disinfected after each use. Alcohol-based surface disinfecting agents are to be used. Please heed hereby the instructions and reaction times given by the manufacturer.

**Maintenance:** Damaged parts must be removed and/or replaced.

**Check and performance test:** Check moving parts such as the dampers for smooth operation. Visually inspect instrument for damage and wear. Before use, ensure that all screws are tightened.

**Packaging:** packaged individually: a standard polyethylene bag can be used. The size of the bag must correspond to the size of the instrument so that it closes properly without straining the bag closure.

**Sterilisation:** Not a surgical instrument. No sterilisation required/ possible.

**Storage:** No special requirements.



**Disposal:** The product can be disposed of at a location that accepts scrap metal.

**Common accessories:** Further information about these articles can be found on our homepage: [www.kawemed.de](http://www.kawemed.de).

Manufacturer: **KaWe**

**Contacting the manufacturer:** Address or telephone number of your dealer or dial +49-7141-68188-0.

**Symbol key:**

	Manufacturer
	Date of manufacture
	Lot code
	Heed the User's Manual
	Caution!
	Product reference number
	Meets EC Standards for medical products

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**KIRCHNER & WILHELM** GmbH + Co. KG