
BODY COMPOSITION ANALYZER

SC-240 / SC-240IM

Instruction manual

<Usage Conditions>

Temperature Range for Use	: 0°C – 35°C (32°F – 95°F)
Relative Humidity	: 30% – 80% (without condensation)
Max Altitude	: 2,000m ASL (6,500ft ASL)
Pressure Range of Environment	: 86kPa – 106kPa

<Storage Conditions>

Temperature Range of Environment	: -10°C – 60°C (14°F – 140°F)
Range of Relative Humidity	: 10% – 90% (without condensation)
To avoid malfunctions, avoid storing the equipment where there is direct sunlight, significant temperature changes, the risk of dampness, a large amount of dust, in the vicinity of fires, or where there is the risk of receiving vibrations or shocks.	

Please read this Instruction Manual carefully and keep it for future reference.

Intended use

- This equipment can be used in the screening of certain adult diseases and conditions related to body weight and composition.
- It can be used in the monitoring and prevention of conditions caused by excessive deposits of fat tissue such as diabetes, hyperlipidemia, cholelithiasis and fatty liver.
- It can be used in the monitoring of changes in individuals' body composition related to differences in the ratio of fat tissue to lean.
- It can be used to assess the effectiveness of individuals' nutrition and exercise programmers, both for health and physical fitness.
- The TANITA Body Composition Analyzer is indicated for use in the measurement of weight and impedance , and the estimation of body mass index (BMI), total body fat percent, total body water percent and weight , muscle mass (skeletal and smooth), physique rating, bone mass, visceral fat rating with healthy range, basal metabolic rate (BMR), and metabolic age, using B IA (Bioelectrical Impedance Analysis).

The device is indicated for use for healthy children 5-17 years old and healthy adults with active, moderately active, to inactive lifestyles.

Efficacy

1. This product is simple to use, and requires no specialized facilities or expertise to take measurements.
2. Measurements can be taken quickly and easily, causing minimal inconvenience to the patient during measurements.

Contents


Before use (safety precautions)	Safety Notes	3
	Part Names & Accessories	6
	Names and Functions of Display Panel & Operation Keys	7
	Preparation (Power supply)	8
	Various setting	9
How to use (cautions for safety)	How to use (Mode selection)	10
	How to use (Body composition mode)	11
	How to use (Scale mode)	14
	How to use (Scale / Weight lock mode)	15
	How to use (Scale / BMI mode)	16
When necessary (cautions for safety)	Troubleshooting	17
	Various criteria	18
	Connection with a personal computer	22
	Technical notes	26
	Specifications	30

How to Measure Correctly

The following outlines precautionary measures to be taken to avoid injury to the users of this device and others, and to prevent damage to property. Please familiarise yourself with the contents to ensure the safe operation of this equipment.

	Warning Ignoring instructions highlighted with this mark could result in death or severe injury.
	Caution Ignoring instructions highlighted with this mark could result in injury or damage to property.
	This mark indicates actions that are prohibited.
	This mark indicates instructions that must always be adhered to.

Warning

 Prohibited	<p>This device must not be used on subjects with pacemakers or other mechanical implants.</p> <p>This device passes a weak electrical current through the body which could interfere with and cause the malfunction of electrical medical implants, resulting in serious repercussions.</p>	
	<p>Do not handle the plug with a wet hand.</p> <p>Such action may result in electric shock, fire, or leakage.</p>	
	<p>Keep this equipment away from flammable gas and oxygen rich environment.</p>	
	<p>No modification of this equipment is allowed.</p> <p>There is a risk of electrocution or injury. Nor can precise analysis be guaranteed.</p>	
	<p>Do not use double adapters.</p> <p>May result in fire.</p>	
 Always...	<p>Use only a correctly wired outlet.</p> <p>Only genuine cables and equipments are allowed.</p>	

Caution

Keep off wet circumstance.

Avoid using on subjects with allergies to metals.

Allergic reactions may be caused by the stainless steel used in the electrodes of this device.

Do not jump on the equipment.

Do not lean on the equipment.

Do not use this equipment near other products that emit electromagnetic waves.

Do not insert fingers into gaps and holes.

Avoid impacts on the display.

The screen panel may split and cause injury.

Ensure disabled persons do not attempt to use the equipment by themselves.

A companion must perform the measurements.

Prohibited

Be sure to clean the scale platform with appropriate disinfectant after each use.

For your safety to avoid the risk of electric shock and to secure the accuracy, keep a clearance with patients during measuring.

Continually monitor both the subject and the equipment for anomalies.

Should an anomaly with the subject or the equipment be discovered, take appropriate action, such as stopping the equipment, while ensuring the safety of the subject.

Unplug the equipment from the wall outlet when it will not be in use for long periods of time.

Do not interpret analysis results yourself (including evaluating measurements and formulating exercise programmes based on results, etc.).

Weight loss and exercise based on self-analysis could be detrimental to you health. Always follow the advice of a qualified professional.

This equipment is designated a Class B IT equipment (mainly for systems intended to be used in domestic environments) and CE (EMC), but weak equipment may be affected.

If you connect your computer or peripherals, please use the compliance with IEC60601-1 (EN60601-1). Or if using the compliance with IEC60950 (EN60950), power must be supplied from themedical isolating transformer. Keep at a distance of 1.5m between each equipment during operation. If you ignore this requirement, may cause electrical shock of patients or operation.

Always...

How to Measure Correctly

For Accurate Measurements

Avoid measuring after hard exercise.

May result in measurement errors. Please take measurements after sufficient rest.

Prohibited

Avoid measuring after over-eating or over-drinking, and when severely dehydrated.

May result in measurement errors. For greater accuracy, avoid using directly after waking up. Use at the same time on each occasion, at least three hours after eating.

Do not take measurements while using transmitters, such as mobile phones, which may affect readings.

Use the device under the same conditions and in the same position as much as possible to track changes.

Readings are greatly affected by the level of hydration and position of the body. Please use at the same time each day, under the same conditions and in the same body position.

Always...

Avoid measuring in locations of greatly differing temperature.

May result in measurement errors. When the equipment has been moved to a location $\pm 20^{\circ}\text{C}$ ($\pm 36^{\circ}\text{F}$), allow to stand for at least 2 hours before using.

Bare feet should be placed correctly on the electrode panel for measurement.

Also, make sure the soles of feet are free of excess dirt, as this may also act as a barrier to the mild current.

Use in a stable location.

Errors in measurements may occur when the device is used in an unstable location.

Scheduled Maintenance

Perform periodic checks TANITA recommends that each facility conduct periodic checks of each unit.

1. Check the following at least daily:
 - Place the unit on a stable and level surface.
2. Visually inspect the following at least weekly:
 - the display for any damage or contamination
 - all cables, cords, and connector ends for damage or contamination
 - all safety-related labeling for legibility
 - all accessories (sensors, electrodes) for wear or damage
 - documentation for current revision of the unit
3. Visually inspect the following at least monthly:
 - Screws on unit


Always...

Update settings, replace items, or call for service as necessary based on results of visual inspection. Do not use the equipment if you see any signs of damage. Qualified service personnel must check any equipment that is damaged for proper operation before putting the equipment back into use.

The equipment should only be inspected in accordance with the regulations in each country.

Names and Functions of Display Panel & Operation Keys



- ON/OFF** ON / OFF key : turn the power ON / OFF
- Mode** Mode selection key : Select Body composition analyzer / Weigh scale
- PT/CE** PT / CE : Set preset tare / Clear input values
- ▲** UP : Increase numerical values
- ▼** DOWN : Decrease numerical values
- Enter/Zero** Enter / Zero reset : Enter input value / Scale reset to zero / Mode setting
- NET** : Indicates that the clothes weight is input
- Weight Only** : Indicates that the Weight Only (Scale) Mode
- Weight Lock** : Indicates that the Weight is locking
- Weight Lock Mode** : Indicates that the Weight Lock function is activated
- Age** : Indicates that the age is input
- Height** : Indicates that the height is input
- Step On** : Indicates that the measuring start

<Symbols and their Meanings>

P T	Preset tare (Clothes weight)		Male		Female
	Athletic mode		Direct current		Input, Output
	Caution. Refer to the attached information		Alternating current		Negative polarity
	⊕ for positive polarity				

Preparation (Power supply)


<Using batteries>

- Please change the batteries LR6 (AA) carefully to avoid dropping them on your feet.
- Ensure that the batteries are inserted with the correct polarity \oplus/\ominus . If the polarity is incorrect, then the batteries may leak and damage the equipment.
- When not in use for a long time, remove the batteries before storing the equipment.

<Using the AC adapter>

- (1) Insert the AC adapter jack into the AC adapter inlet on the right side of the main unit.
- (2) Plug the AC adapter into the power outlet.

⚠ Caution

- This equipment must be used with the included AC adapter.
- Please observe the following instructions for accurate measurement.

Measurement may not be possible on the unstable environment.


During weight measurement, please don't touch any connecting cable "such as a AC Adapter cable and PC communication cable" to avoid causing unstable scale installation.

Various setting

Call up the setting mode.

1 Press ON/OFF to turn on the power.

2 Press Enter/Zero for 1 second.

SET 0 is displayed.

3 Select the setting items.

- 1) by ▲ or ▼,
- 2) and Enter/Zero key.

SET1 Set ON / OFF of the beep sound.

(0.off or 1.on, default: 1.on)

SET2 Set ON / OFF of the athlete selection.

(0.off or 1.on, default: 1.on)

SET3 Set the automatic acceptance time when inputting.

(input range 0-9 seconds, default: 5 seconds)

* "0" automatic acceptance function deactivate.


SET4 Set the automatic power off time.

(0, 5, 10, 30, 60 minutes, default: 5 minutes)

* "0" automatic power off function deactivate.


SET5 Set the unit.

(0kg or 1.lb, default: 1.lb)

4 Input the setting value.


- 1) by ▲ or ▼,
- 2) and Enter/Zero key.

5 Press ON/OFF to turn the power off.

How to use (Mode selection)

1 Press ON/OFF to turn on the power.


Model name is displayed.


00 is displayed (Body composition analyzer).

Mode is selected by **Mode** key.


Body composition mode (P.11)


Scale mode (P.14)

Scale / BMI mode (P.16)

Scale / Weight lock mode (P.15)

How to use (Body composition mode)

1 Press ON/OFF to turn on the power.



2 Check that the Body composition analyzer mode is selected (☞ P.10).

3 Input the clothes weight.


1) by ▲ or ▼,

2) and Enter/Zero key


* PT/CE ➡ clear the input.

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

4 Select the Body type and Gender.


1) by ▲ or ▼,

2) and Enter/Zero key


* PT/CE ➡ return to previous status.

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero



How to use (Body composition mode) (continued)

5 Input the age

1) by ▲ or ▼,

2) and **Enter/Zero** key

* **PT/CE** → return to previous status



6 Input the height

1) by ▲ or ▼,






2) and **Enter/Zero** key

* **PT/CE** → return to previous status

How to use (Body composition mode)

7 Step on the platform with bare feet, after “Step on” flashing.

8 Measurement completion.

* PT/CE, Enter/Zero → Clear the displayed result

How to use (Scale mode)

1 Press ON/OFF to turn on the power.





2 Check that the Weight only mode is selected (☞ P.10).

3 Input the clothes weight.

1) by ▲ or ▼,



2) and Enter/Zero key

* PT/CE → clear the input.

4 Step on the platform, after “Step on” flashing.

5 Measurement completion.

How to use (Scale / Weight lock mode)

Activation of Weight lock function.

- 1 Press ON/OFF to turn on the power.
- 2 Check that the Weight only mode (Weight lock mode) is selected (P.10).

3 Input the clothes weight.

1) by ▲ or ▼,



2) and Enter/Zero key

* PT/CE → clear the input.


4 Step on the platform, after “Step on” flashing.

5 Measurement completion.



* Enter/Zero → Clear the displayed result

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

How to use (Scale / BMI mode)

1 Press ON/OFF to turn on the power.



2 Check that the BMI mode is selected
( P.10).

3 Input the clothes weight.



1) by ▲ or ▼,

2) and Enter/Zero key


* PT/CE → clear the input.

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

4 Input the height



1) by ▲ or ▼,

2) and Enter/Zero key


* PT/CE → return to previous status

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero

ON/OFF Mode PT/CE ▲ ▼ Enter/Zero


5 Step on the platform, after “Step on” flashing.

6 Measurement completion.

* PT/CE, Enter/Zero → Clear the displayed result

- Please check the following before asking for repair.

	Symptom	Please check
How to measure	<p>Impedance measurement error</p> 	<ul style="list-style-type: none"> • Measure with bare feet. • When the soles of your feet are dry, drop water for about 0.5 ml with the attached dropper on the electrodes before measurement. • Check the input information.
	<p>Zero point error</p> 	<ul style="list-style-type: none"> • Turn off the power, and remove the items on the platform and turn on the power again, and then redo the measurement.
	<p>The measured weight is not stable.</p>	<ul style="list-style-type: none"> • Is it installed at a place with vibrations? • Is the platform inclined? ⇒ Keep the platform horizontal. • Is anything caught in the gaps of the platform? ⇒ Remove anything caught in the gaps.
Display part	<p>Nothing is displayed even after turning on the power.</p>	<ul style="list-style-type: none"> • Check that the power supply is connected correctly.
	 <p>is displayed.</p>	<ul style="list-style-type: none"> • The weight to measure exceeds the measurement range.

Note

• In the case that an error is displayed other than the above, turn off the power once, and then measure again.

If the same error is displayed repeatedly, contact our customer service center.

- What is total body water percentage? (Applicable age 18 - 99)

Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight.

Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage will ensure the body functions efficiently and will reduce the risk of developing associated health problems.

Your body water levels naturally fluctuate throughout the day and night. Your body tends to be dehydrated after a long night and there are differences in fluid distribution between day and night. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may cause variations in your hydration levels.

Your body water percentage reading should act as a guide and should not be used to specifically determine your absolute recommended total body water percentage. It is important to look for long-term changes in total body water percentage and maintain a consistent, healthy total body water percentage.

Drinking a large quantity of water in one sitting will not instantly change your water level. In fact, it will increase your body fat reading due to the additional weight gain. Please monitor all readings over time to track the relative change.

Every individual varies but as a guide the average total body water percentage ranges for a healthy adult are:

Female : 45 to 60%

Male : 50 to 65%

Source : Based on Tanita's Internal Research

Note: The total body water percentage will tend to decrease as the percentage of body fat increases. A person with a high percentage of body fat may fall below the average body water percentage. As you lose body fat the total body water percentage should gradually move towards the typical range given above.

- What is visceral fat rating? (Applicable age 18 - 99)

This feature indicates the rating of visceral fat in your body.

Visceral fat is the fat that is in the internal abdominal cavity, surrounding the vital organs in the trunk (abdominal) area. Research shows that even if your weight and body fat remains constant, as you get older the distribution of fat changes and is more likely to shift to the trunk area especially post menopause. Ensuring you have healthy levels of visceral fat may reduce the risk of certain diseases such as heart disease, high blood pressure, and the onset of type 2 diabetes.

The Tanita Body Composition Analyzer will provide you with a visceral fat rating from 1 – 59.

Rating from 1 to 12

Indicates you have a healthy level of visceral fat. Continue monitoring your rating to ensure that it stays within the healthy range.

Rating from 13 – 59

Indicates you have an excess level of visceral fat. Consider making changes in your lifestyle possibly through diet changes or increasing exercise.

Source: 1) Tanita Institute Contract Study. Algorithm Development for Estimating Visceral Fat Rating. SB Heymsfield MD. Columbia University College of Physicians and Surgeons 2004.

- 2) Wang, Z., et al. Japanese-American Differences in Visceral Adiposity and a Simplified Estimation Method for Visceral Adipose Tissue. North American Association for the Study of Obesity. Annual Meeting. Abstract 518-P. 2004

Note:

- Even if you have a low body fat rate, you may have a high visceral fat level.
- For medical diagnosis, consult a physician.

Various criteria (continued)

- What is basal metabolic rate (BMR)? (Applicable age 18 - 99)

WHAT IS BMR?

Your Basal Metabolic Rate (BMR) is the minimum level of energy your body needs when at rest to function effectively including your respiratory and circulatory organs, neural system, liver, kidneys, and other organs. You burn calories when sleeping.

About 70% of calories consumed every day are used for your basal metabolism. In addition, energy is used when doing any kind of activity however; the more vigorous the activity is the more calories are burned. This is because skeletal muscle (which accounts for approximately 40% of your body weight) acts as your metabolic engine and uses a large amount of energy. Your basal metabolism is greatly affected by the quantity of muscles you have, therefore increasing your muscle mass will help increase your basal metabolism.

By studying healthy individuals, scientists have found that as people age, their metabolic rate changes. Basal metabolism rises as a child matures. After a peak at the age of 16 or 17, it typically starts to decrease gradually.

Having a higher basal metabolism will increase the number of calories used and help to decrease the amount of body fat. A low basal metabolic rate will make it harder to lose body fat and overall weight.

HOW DOES A TANITA BODY COMPOSITION ANALYZER CALCULATE BMR?

The basic way of calculating Basal Metabolic Rate BMR is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).*

* Reliability on equation for Basal Metabolic Rate: At 2002 Nutrition Week : A Scientific and Clinical Forum and Exposition Title: International Comparison: Resting Energy Expenditure Prediction Models: The American Journal of Clinical Nutrition

- What is metabolic age? (Applicable age 18 - 99)

This feature calculates your BMR and indicates the average age associated with that type of metabolism.

If your BMR Age is higher than your actual age, it is an indication that you need to improve your metabolic rate. Increased exercise will build healthy muscle tissue, which will improve your metabolic age.

You will obtain a reading between 12 and 90. Under 12 will be displayed as "12" and over 90 displayed as "90".

- What is muscle mass? (Applicable age 18 - 99)

This feature indicates the weight of muscle in your body. the muscle mass displayed includes the skeletal muscles, smooth muscles (such as cardiac and digestive muscles) and the water contained in these muscles.

Muscles play an important role as they act as an engine in consuming energy. As your muscle mass increase, your energy consumption increases helping you reduce excess body fat levels and lose weight in a healthy way.

- What is bone mass? (Applicable age 18 - 99)

This feature indicates the amount of bone (bone mineral level, calcium or other minerals) in the body.

Research has shown that exercise and the development of muscle tissue are related to stronger, healthier bones. While bone structure is unlikely to make noticeable changes in a short period, it is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercise. People worried about bone disease should consult their physician. People who suffer from osteoporosis or low bone densities due to advanced age, young age, pregnancy, hormonal treatment or other causes, may not get accurate estimations of their bone mass.

Below is the result of estimated bone masses of persons aged 20 to 40, who are said to have the largest amounts of bone masses, by weight. (Source : Tanita Body Weight Science Institute)

Please use the below charts as a guide to compare your bone mass reading.

Women : Average of estimated bone mass

Weight (lb)		
Less than 110 lb	110lb - 165 lb	165 lb and up
4.3 lb	5.3 lb	6.5 lb

Weight (kg)		
Less than 50 kg	50 kg - 75 kg	75 kg and up
1,95 kg	2,40 kg	2,95 kg

Men : Average of estimated bone mass

Weight (lb)		
Less than 143 lb	143 lb - 209 lb	209 lb and up
5.9 lb	7.3 lb	8.1 lb

Weight (kg)		
Less than 65 kg	65 kg - 95 kg	95 kg and up
2,66 kg	3,29 kg	3,69 kg

Note: Persons described below may obtain varying readings and should take the values given for reference purposes only.

- Elderly persons
- Women during or after menopause
- People receiving hormone therapy

“Estimated bone mass” is a value estimated statistically based on its correlation with the fat-free amount (tissues other than the fat). “Estimated bone mass” does not give a direct judgment on the hardness or strength of the bones or the risks of bone fractures. If you have concerns over your bones, you are recommended to consult a specialist physician.

Connection with a personal computer

The USB interface enables input and output from this equipment.

- Specifications

Communication standards	USB
Communication method	Asynchronous Transmission
Signal speed	9600 baud
Data bit length	8 bits
Parity	None
Stop bit	1 bit
Flow control	None
Terminator	CR+LF

Note

- USB connector (B - Type 4 - pin female) are located on the back of the control box.
- Please provide your own cable as necessary as none are included.
USB cable : A - Type 4 - pin (male) - B - Type 4 - pin (male)
- You must install the necessary driver onto your PC, available download from <http://www.tanita.com> or on the enclosed CD-Rom.

- Transmit data

Caution

Transmission Data is output immediately after measurement regardless of the status of the receiving equipment (personal computer, etc.). Therefore the receiving equipment needs to be ready to accept the data before measuring.

(1) Output example (in the case to kg unit)

Note

- The data are divided with commas (,) for each data.
- The terminator (end of the data) is CR (ASCII code 0DH), LF (ASCII code 0AH).
- The control data for the items 1 – 4 are for expansion. They are not used presently so the receiving side can ignore them.
- Because of its specifications, the scale will output an irrelevant signal approximately 0.1 second after data output. Make sure you ignore this signal and do not import it.

Entire body data

Model number		Serial number		Body type		Gender		Age	
MO	"XXXXXX"	SN	"XXXXXXXX"	Bt	0or2	GE	1or2	AG	XX
Height		Preset tare		Body weight		Fat %		Fat mass	
Hm	XXX.X	Pt	XX.X	Wk	XXX.X	FW	XX.X	fW	XXX.X
Fat free mass		Muscle mass		Muscle score		Bone mass		Body water	
MW	XXX.X	mW	XXX.X	sW	XX	bW	XXX.X	wW	XXX.X
Body water %		BMI		Ideal body weight		Degree of obesity		Visceral fat rating	
ww	XX.X	MI	XXX.X	Sw	XXX.X	OV	XXX.X	IF	XX
BMR (kJ)		BMR (kcal)		BMR score		Metabolic age		Impedance (50kHz)	
rb	XXXXX	rB	XXXXX	rJ	XX	rA	XX	ZF	XXXX.X
Checksum									
CS	XX								

Connection with a personal computer (continued)

(2) Output data items

<kg / cm mode>

Item	Header	format	kg / cm mode	Order of output				
				Body composition			Scale	BMI
				Adult	Athlete	Child		
Control data	{0		2 bytes fixed length	1	1	1	1	1
Control data	~0		1 byte fixed length	2	2	2	2	2
Control data	~1		1 byte fixed length	3	3	3	3	3
Control data	~2		1 byte fixed length	4	4	4	4	4
Model number	MO	"XXXXXX"	8 bytes fixed length ("SC-240")	5	5	5	5	5
Serial number	SN	"XXXXXXXX"	10 bytes fixed length	6	6	6	6	6
Body type	Bt	0 or 2	1 byte fixed length (0: Standard 2: Athletic)	7	7	7		
Gender	GE	1 or 2	1 byte fixed length (1: Male 2: Female)	8	8	8		
Age	AG	XX	1-2 bytes variable length	9	9	9		
Height	Hm	XXX.X	4-5 bytes variable length (cm)	10	10	10		7
Preset tare	Pt	XX.X	3-4 bytes variable length (increments 0.1 kg)	11	11	11	7	8
Body weight	Wk	XXX.X	3-5 bytes variable length (increments 0.1 kg)	12	12	12	8	9
Fat %	FW	XX.X	3-4 bytes variable length (increments 0.1%)	13	13	13		
Fat mass	fW	XXX.X	3-5 bytes variable length (increments 0.1 kg)	14	14	14		
FFM	MW	XXX.X	3-5 bytes variable length (increments 0.1 kg)	15	15	15		
Muscle mass	mW	XXX.X	3-5 bytes variable length (increments 0.1 kg)	16	16			
Muscle score	sW	XX	1-2 bytes variable length (1-24)	17	17			
Bone mass	bW	XXX.X	3-5 bytes variable length (increments 0.1 kg)	18	18			
Body water	wW	XXX.X	3-5 bytes variable length (increments 0.1 kg)	19	19			
Body water %	ww	XX.X	4 bytes fixed length (increments 0.1%)	20	20			
BMI	MI	XXX.X	3-5 bytes variable length (increments 0.1)	21	21	16		10
Ideal body weight	Sw	XXX.X	3-5 bytes variable length (increments 0.1 kg)	22				
Degree of obesity	OV	XXX.X	3-5 bytes variable length (increments 0.1%)	23				
Visceral fat rating	IF	XX	1-2 bytes variable length	24	22			
BMR (kJ)	rb	XXXXX	1-5 bytes variable length (increments 1 kJ)	25	23			
BMR (kcal)	rB	XXXXX	1-5 bytes variable length (increments 1 kcal)	26	24			
BMR score	rJ	XX	1-2 bytes variable length (1-24)	27	25			
Metabolic age	rA	XX	2 bytes fixed length	28	26			
Impedance (50 kHz)	ZF	XXXX.X	1-2 bytes variable length (4-55)	29	27	17	9	
Checksum	CS	XX	2 bytes fixed length	30	28	18	10	11

Output sample

Standard, Male, 30 years old, Height 170 cm

{0, 16, ~0, 1, ~1, 1, ~2, 1, MO, "SC-240", SN, "00000002", Bt, 0, GE, 1, AG, 30, Hm, 170, Pt, 0.0, Wk, 82.7, FW, 18.6, fW, 15.4, MW, 67.3, mW, 64.0, sW, 19, bW, 3.3, wW, 48.0, ww, 58.0, MI, 28.6, Sw, 63.6, OV, 30.0, IF, 6, rb, 8255, rB, 1973, rJ, 12, rA, 28, ZF, 369.4, CS, FD

<lb / ft-in mode>

Item	Header	format	lb / ft-in mode	Order of output				
				Body composition			Scale	BMI
				Adult	Athlete	Child		
Control data	{0		2 bytes fixed length	1	1	1	1	1
Control data	~0		1 byte fixed length	2	2	2	2	2
Control data	~1		1 byte fixed length	3	3	3	3	3
Control data	~2		1 byte fixed length	4	4	4	4	4
Model number	MO	"XXXXXX"	8 bytes fixed length ("SC-240")	5	5	5	5	5
Serial number	SN	"XXXXXXXX"	10 bytes fixed length	6	6	6	6	6
Body type	Bt	0 or 2	1 byte fixed length (0: Standard 2: Athletic)	7	7	7		
Gender	GE	1 or 2	1 byte fixed length (1: Male 2: Female)	8	8	8		
Age	AG	XX	1-2 bytes variable length	9	9	9		
Height	Hm	XXX.X	4 bytes variable length (in)	10	10	10		7
Preset tare	Pt	XX.X	3-4 bytes variable length (increments 0.2lb)	11	11	11	7	8
Body weight	Wk	XXX.X	3-5 bytes variable length (increments 0.2 lb)	12	12	12	8	9
Fat %	FW	XX.X	3-4 bytes variable length (increments 0.1%)	13	13	13		
Fat mass	fW	XXX.X	3-5 bytes variable length (increments 0.2 lb)	14	14	14		
FFM	MW	XXX.X	3-5 bytes variable length (increments 0.2 lb)	15	15	15		
Muscle mass	mW	XXX.X	3-5 bytes variable length (increments 0.2 lb)	16	16			
Muscle score	sW	XX	1-2 bytes variable length (1-24)	17	17			
Bone mass	bW	XXX.X	3-5 bytes variable length (increments 0.2 lb)	18	18			
Body water	wW	XXX.X	3-5 bytes variable length (increments 0.2 lb)	19	19			
Body water %	ww	XX.X	4 bytes fixed length (increments 0.1%)	20	20			
BMI	MI	XXX.X	3-5 bytes variable length (increments 0.1)	21	21	16		10
Ideal body weight	Sw	XXX.X	3-5 bytes variable length (increments 0.2 lb)	22				
Degree of obesity	OV	XXX.X	3-5 bytes variable length (increments 0.1%)	23				
Visceral fat rating	IF	XX	1-2 bytes variable length	24	22			
BMR (kJ)	rb	XXXXX	1-5 bytes variable length (increments 1 kJ)	25	23			
BMR (kcal)	rB	XXXXX	1-5 bytes variable length (increments 1 kcal)	26	24			
BMR score	rJ	XX	1-2 bytes variable length (1-24)	27	25			
Metabolic age	rA	XX	2 bytes fixed length	28	26			
Impedance (50 kHz)	ZF	XXXX.X	5-6 bytes variable length (increments 0.1)	29	27	17	9	
Checksum	CS	XX	2 bytes fixed length	30	28	18	10	11

Output sample

Standard, Male, 30 years old, Height 67.0 in

{0, 16, ~0, 1, ~1, 1, ~2, 1, MO, "SC-240", SN, "00000003", Bt, 0, GE, 1, AG, 30, Hi, 67.0, Pt, 0.0, Wp, 182.2, FW, 18.7, fW, 34.0, MW, 148.2, mW, 140.8, sW, 19, bW, 7.4, wW, 105.6, ww, 58.0, MI, 28.5, Sw, 140.4, OV, 29.8, IF, 6, rb, 8247, rB, 1971 ,rJ, 12, rA, 29, ZF, 372.8, CS, 2A

Technical notes

Body composition measurement by the BIA method.

Introduction

This equipment provides estimated values for each measured value of body fat percentage, fat mass, fatfree mass, muscle mass and bone mass by the DXA method, estimated value for the total body water measured value by the dilution method and estimated value for the visceral fat rating by MRI method using the Bioelectrical Impedance Analysis (BIA method).

For measurement, a mode must be selected based on body type.

- 1) Standard (for 5-99 years of age)
- 2) Athletic (for Athletic persons who exercise considerably more than non-athlete)

Making a distinction by body type in the measurement mode produces more reliable body composition measurements for athletic persons, whose body compositions differ from those of average persons.


- Principles of body composition measurement

BIA is a means of measuring body composition – fat mass, predicted muscle mass, etc. – by measuring bioelectrical impedance in the body. Fat within the body allows almost no electricity to pass through, while electricity passes rather easily through water, much of which is found in muscles. The degree of difficulty with which electricity passes through a substance is known as the electrical resistance, and the percentage of fat and other body constituents can be inferred from measurements of this resistance.

The Tanita Body Composition Analyzer measures body composition using a constant current source with a high frequency current (50 kHz, 90 μ A). The 8 electrodes are positioned so that electric current is supplied from the electrodes on the tips of the toes of both feet, and voltage is measured on the heel of both feet. The current flows into the upper limbs or lower limbs, depending on the body part(s) to be measured.

- What is the DXA method?

DXA was originally designed to measure bone mineral content, but in the full-body scan mode the body fat percentage, fat mass, and fat free mass of individual body parts (arms, legs, trunk) can also be measured. The image below shows one example of body composition measurement results obtained by DXA.

Body composition measurement results obtained by DXA (Lunar Co., Ltd; DPX-L)

- What is dilution method?

In the dilution method, a labeled substance for a known amount is given and the concentration in equilibrium diffusing evenly is measured to obtain the total amount of the solvent that dilutes the labeled substance.

To measure the total body water (TBW), deuterium oxide (D₂O) is generally used as the labeled substance. Deuterium oxide uses the overall total body water as dilution space so the total body water can be obtained. To obtain the extracellular fluid amount, sodium bromide (NaBr) is used as a labeled substance. Bromine (Br) is said to not enter the inside of cells, and uses extracellular fluid as the dilution space.

- What is the visceral fat?


Visceral Adipose Tissue (VAT) has been associated with increased risk of developing lifestyle-related diseases. Accordingly, knowing and periodically checking the estimated VAT accumulation serves as one factor of a number of factors in assessing the prevention of lifestyle related diseases.

Tanita has developed the technology for measuring the VAT accumulation risk through bioelectrical impedance analysis (BIA) in comparison with image analysis applied to magnetic resonance imaging (MRI), in addition to the established technology for measuring the percent of body fat. The VAT accumulation risk is calculated by estimating the VAT area by the BIA method on the basis of MRI image processing. This method has a higher correlation than the estimation of the VAT accumulation risk based on BMI or abdominal circumference (waist circumference), allowing estimation that corresponds more precisely to individuals.


* The VAT area by MRI is calculated by carrying out an image processing of the cross section of the lumbar vertebra L4-L5 regions.

(Fig. 1 - Fig. 3: Research results by N. Y. Columbia University and Jikei University Published by the North American Association for the Study of Obesity [NAASO] in 2004.)


<Fig. 1> Relationship between VAT Area and BMI

<Fig. 2> Relationship between VAT Area and Waist Circumference

<Fig. 3> Relationship between VAT Area by MRI and Estimated VAT Area by Tanita's BIA

Technical notes (continued)

- Factors giving errors in measurement

In the BIA method, impedance is measured and the body composition is calculated based on the value. It is known that impedance changes by the amount of the total body water that occupies about 60% of weight and the change in its distribution and temperature change. Therefore, for the purpose of research or for daily repeating of measurements, the measurement conditions must be kept constant. Measurement under the changing conditions of temperature and total body water distribution or blood flow volume of extremities due to exercising, taking a bath, etc., affects the measurement result since the electric resistance in the body also changes.

Therefore, it is recommended to measure under the following conditions for stable measurement.

- 1) 3 hours have passed after getting up and normal lifestyle activities are carried out during this period. (The impedance transits staying at a high level if you remain sitting after getting up or drive a car, etc.)
- 2) 3 hours or more have passed after eating. (For 2 – 3 hours after eating, the impedance has a tendency to decrease.)
- 3) 12 hours or more have passed after vigorous exercise for measurement. (The tendency toward changes in impedance is not stable depending on the type and rigorousness of the exercise.)
- 4) If possible urinate before taking measurement.
- 5) For repeated measurements, measure at the same hour as much as possible. (At the same time of measurement of weight, the measurements can be made more stable by measuring at the same time of the day)

Very stable measured values can be obtained by measuring under the above conditions.

And in the development of this equipment, the following 6 items were set as conditions for the regression equation.

- 1) Prohibition of alcohol intake for 12 hours before measurement
- 2) Prohibition of excessive exercise for 12 hours before measurement.
- 3) Prohibition of excessive eating and drinking the day before measurement
- 4) Prohibition of eating and drinking for 3 hours before measurement
- 5) Urination just before measurement
- 6) Avoid measurements during menstruation (for female)

Inter-day changes

The diagrams below offer examples of actual measurements made of inter-day changes. A study was done to determine the degree of change in the impedance between the feet during dehydration; the first two days represent a normal daily routine, while in the latter two days a state of dehydration was induced using a sauna.


No significant inter-day change was measured in body weight, impedance between the feet, or body fat percentage during the normal daily routine. During the dehydrated state, however, a drop in body weight of 1kg was noted, with the impedance between the feet rising approximately 15Ω on the first day of dehydration and $30\text{-}35\Omega$ on the second day. As a result, body fat percentage was up by around 1% on the first day of dehydration and by 1.5-2% on the second day.

As mentioned earlier, impedance increases when body weight is reduced (such as by dehydration), and decreases when body weight is increased through excess consumption of food and drink. The inter-day change in impedance is thus inversely proportional to the change in body weight.

These inter-day changes stem from such causes as:

- 1) Temporary increases in body weight (total body water) through overeating and overdrinking
- 2) Dehydration due to heavy sweating during vigorous exercise
- 3) Dehydration due to alcohol consumption or the use of diuretics
- 4) Dehydration due to heavy sweating during saunas, etc.

Accordingly, it is recommended that instructions be provided to the subject to help eliminate these causes when accurate measurements are needed.

Specifications

Model		SC-240 / SC-240IM	
Power source	AC adapter	Centre minus	
	Battery	LR6 (AA) x 6	
Power Consumption		0.5 W	
Impedance Measurement	Measurement System	Tetra polar Bioelectrical Impedance Analysis	
	Measurement Frequency	50 kHz	
	Measurement Current	90 µA	
	Measurement Range	150-1200 Ω	
Weight Measurement	Measurement System	Strain Gauge Load Cell	
	Maximum Capacity	440 lb	200 kg
	Minimum Graduation	0.2 lb	0.1 kg
Input Items	Clothes Weight	0-44 lb / 0.2 lb increments	0-10 kg / 0.1 kg increments
	Gender	Female / Male	
	Body Type	Standard / Athletic	
	Age	5-99 years (Standard) / 18-99 years (Athletic)	
	Height	3 ft-7 ft 11.5 in / 0.5 in increments	90-249 cm / 1 cm increments
Output Items (Display)	Weight	0-440 lb/0.2 lb increments	0-200.0 kg/0.1 kg increments
	Body Fat %	3-75% / 0.1% increments	
	BMI	0.1 increments	
	Total body water %	15 - 85% / 0.1% increments	
Interface Connections		USB	
Weight of Equipment		10.0 lb	4.5 kg
Size		13.4×17.2×2.1 in	341×437×54 mm

USA and Canada

Federal Communications Commission and Canadian ICES Notice

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and Canadian ICES-003. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio or television technician for help.

Modifications

The FCC requires the user to be notified that any changes or modifications made to this device that are not expressly approved by Tanita Corporation may void the user's authority to operate the equipment.

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Visceral adipose tissue (VAT) is fat that accumulates in the abdominal cavity and around internal organs. VAT is said to be more likely to cause lifestyle-related diseases than subcutaneous adipose tissue (SCAT). Accordingly, knowing and periodically checking the VAT accumulation risk serves as an important guide in the prevention of lifestyle-related diseases.

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*The VAT area by MRI is calculated by carrying out an image processing of the cross section of the lumbar vertebra L4-L5 regions.

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TANITA Corporation of America, Inc.
TANITA Corporation
TANITA Europe B.V.
TANITA Health Equipment H.K.LTD.
TANITA India Private Limited
TANITA (Shanghai) Trading Co., Ltd.